

**ENABLE YOUR DEVICE
TO GO ONLINE**

DiagConnector

The **ideal way** to transfer data from a device to a patient data management system



universal

cost-effective

easy-to-handle

UPGR@DE NOW

TECHNICAL DATA

- Universal base station for POC devices (e.g. for opTrilyzer®)
- Usable as software or as independent hardware
- Interfaces: USB, RS232, network
- Export of data according to HL7, LDT, GDT
- Web interface with experimental and patient data management systems & configuration*
- Connection via VPN tunnel feasible*
- Sending of experimental data via UMTS, WIFI (pre-installed)
- CE marking in preparation

*if used as independent hardware



ADVANTAGES

- ✓ fast automatic data transmission from the device to the laboratory data management system
- ✓ maintenance-free operation
- ✓ no moving parts inside
- ✓ automatic software updates if requested

CONTACT

Dr. Frauke Adams
☎ +49 (0) 3302 55199-14

✉ f.adams@diagnostiknet-bb.de
www.diagnostiknet-bb.de